

## **Instructions for using the GenomicLand software**

- ***Step 1:*** Download the R software on the link <https://cran.r-project.org/bin/windows/base.old/3.5.1/> (**version 3.5.1**).
  
- ***Step 2:*** Install the R software as instructed in the topic **Installation and other instructions**.
  
- ***Step 3:*** Set environment variables
  - i. In Search, search for and then select: System (Control Panel);
  - ii. Click the **Advanced system settings** link;
  - iii. Click **Environment Variables**;
  - iv. Under **System variables**, select **Path** and click on **edit**;
  - v. Click **New**, and **add or concatenate ( use ; )** the folder address for R to there (for example, C:\Program Files\R\R-3.5.0\bin\x64);
  - vi. Click **OK**. Close all remaining windows by clicking **OK**.
  
- ***Step 4:*** Fill the form ([https://docs.google.com/forms/d/e/1FAIpQLSfEvdJgrY-lFY8mG-bgv3\\_nccr11P2NRafI60Thga382E2yWA/viewform?usp=sf\\_link](https://docs.google.com/forms/d/e/1FAIpQLSfEvdJgrY-lFY8mG-bgv3_nccr11P2NRafI60Thga382E2yWA/viewform?usp=sf_link)).
  
- ***Step 5:*** You will receive in the email that you informed the link to download.
  
- ***Step 6:*** Download the GenomicLand software.
  
- ***Step 7:*** Unzip the file. Open the extracted folder, find the file GenomicLand.exe, and click Sending a shortcut to the Desktop.

After completing these steps, the GenomicLand software will be ready for use.